



A portable unit designed for facilitating communication between a speaker and wearer/s of T-coil hearing aids. The LoopHEAR™ is ideal for one to one applications such as at service counters, railway stations, banks, in vehicles and at small meetings. This lightweight unit includes, within the casing:- an induction loop, a built-in amplifier and an external earpiece .

FEATURES

- ➔ The LoopHEAR™ can be used either as a portable unit or in a fixed position
- ➔ Built-in microphone
- ➔ On/Off button
- ➔ Low battery indicator
- ➔ Indicator of emission magnetic field
- ➔ Adjustable receiving volume of the additional earpiece



TECHNICAL SPECIFICATIONS

- ➔ Battery time :
8 hours when fully charged
Standby when turned on but not in use: 2 days
Standby when turned off: 4 months
- ➔ 9vDC power adaptor (3 m wire)
- ➔ Size :25 x 22 x 7 cm
- ➔ Weight: 710 g
- ➔ Supplied in a lightweight case

ACCESSORIES INCLUDED

- ➔ External microphone (3.5mm jack)
- ➔ Additional earpiece for non-wearers of hearing aids or for wearer of hearing aids who don't have the T position.
- ➔ Rechargeable battery NiMH
- ➔ Power supply 9VDC
- ➔ Official pictograms "crossed ear" for the signage

OPTIONAL ACCESSORIES

- ➔ LHA3, listener in the form of a handset

